
# Pixelated 2D-Sensing Ionization Chamber

### **Features**

- 42 mm diameter sensitive area
- Compact package
- lonization chamber with 120 pixel readout for position and shape monitoring
- Robust thin FR4 electrodes
- HV loopback
- Compatible with I128 and I6400 readout electronics




Applications	<ul> <li>Particle therapy isocenter diagnostic systems</li> <li>Beam shape, position and trajectory monitoring</li> <li>General high energy ion beam diagnostics</li> </ul>
Options	Thin film electrode version available (PX-1)

## **Specifications**

Beam compatibility	
Species	Protons, deuterons, fully-stripped carbon
Energy range	30 MeV/nucleon to 500 MeV / nucleon
Beam current density range	Up to 20 nA cm <sup>-2</sup> (particle current)

Sensor	
Туре	Parallel plate single-gap ionization chamber with pixelated cathode
High voltage	500-1000 V nominal (1660 to 3330 V cm <sup>-1</sup> ); maximum 1500 V
Sensitive area	42 mm nominal diameter




Sensor (cont)	
Sensitive volume	Active volume: Pixelated cathode to anode. 7.62 mm spacing.
Pixel geometry	120 pixels 3.80 mm pitch (50 µm inter-pixel gaps typical)
Gain uniformity	Better than +/-2% for beams within the sensitive area.
Position accuracy	Integral linearity better than 50 µm maximum deviation relative over the sensitive area.
Position resolution	Depends on signal to noise ratio; 10's of $\mu m$ achievable provided beam covers more than one strip.
Fiducials	Electrode pixel position tolerance build-up relative to fiducial features on body +/- $0.3$ mm nominal, < +/- $0.1$ mm typical .

Chamber gas	
Operating gas	Dry atmospheric air


Mechanical				
Insertion length	32 mm window to window, 37 mm housing face to face.			
Overall size	170 mm by 170 mm by 37mm approx (see figures)			
Weight	1.3 kg ( 2.8 lb) excluding any added mounting brackets.			
Operating environment	Clean and dust-free, 0 to 35 C (15 to 25 C recommended , < 70% humidity, non-condensing, vibration < 0.1g all axes (1 to 50 Hz) Temperature and pressure compensation of chamber gain must be performed.			
Shipping and storage environment	-10 to 50 C, < 80% humidity, non-condensing, vibration < 1g all axes, 1 to 20 Hz			

## **Beam scattering** Layers in beam path 1 12.5 μm Polyimide foil aluminized both sides 0.1 µm (window) 2 15.9 mm | Air (non-active gap) 152 μm FR4 with copper patterning 5 µm both sides (cathode) 4 7.62 mm | Air (active gap) 5 12.5 μm | Polyimide foil aluminized both sides 0.1 μm (anode) 6 Air (non-active gap) 8.4mm 12.5 μm Polyimide foil aluminized both sides 0.1 μm (window)



Cathode pixel pattern faces the anode.

## Orientation and pixel pattern




The nominal beam entry face for pixel identification is indicated, although the beam can pass through in other direction also.

In the orientation shown the beam passes through the cathode then the anode. The view of the cathode pixel pattern on the next page is looking in the direction of the arrow. The pattern is on the downstream side of the cathode, facing the anode.

### Orientation and pixel pattern

View along nominal beam direction as shown in the previous figure. View through the cathode to the pixel pattern on the downstream side.

Note: there are 120 pixels, numbered from 0 upwards. Some numbers (30, 31, 62, 63, 94, 95) are absent.



#### Connectors Pixel readout Four DSub male high density 44 pin. HD44 #1. Pixel 28 Test pixel 16 n/c 31 1 2 Pixel 27 17 Pixel 29 32 n/c 3 Pixel 25 18 Pixel 26 33 AGnd / KGnd 4 Pixel 23 19 Pixel 24 34 AGnd / KGnd 5 Pixel 21 20 Pixel 22 35 AGnd / KGnd Pixel 19 Pixel 20 6 21 36 AGnd / KGnd Pixel 17 Pixel 18 7 22 37 AGnd / KGnd Pixel 15 8 23 Pixel 16 38 AGnd / KGnd Pixel 13 Pixel 14 9 24 39 AGnd / KGnd Pixel 11 Pixel 12 AGnd / KGnd 10 25 40 Pixel 9 Pixel 10 AGnd / KGnd 11 26 41 Pixel 7 27 Pixel 8 AGnd / KGnd 12 42 Pixel 5 Pixel 6 n/c 13 28 43 14 Pixel 3 29 Pixel 4 44 Pixel 2 15 Pixel 1 30 Pixel 0 Connector shell is common with ionization chamber body. The pin arrangement is compatible with a pin to pin (M-F) 44-way cable connection to an I6400 or I128 electrometer, with pixel 0 connecting to channel 1 and so on. Test pixels are outside the sensitive area and are available for background noise checks. AGnd is the signal ground designation in the I6400 electrometer. KGnd is the signal ground designation in the I128 electrometer. **CAUTION** Do not expose the device to ionizing radiation beams unless all connections to readout electronics and bias supplies are made, or otherwise grounded. Charge build-up and subsequent arcing damage can occur.

## Connectors (cont)

Pixel readout (cont)

HD44 #2.

1	Pixel 60	16	n/c	31	Test pixel 2
2	Pixel 59	17	Pixel 61	32	n/c
3	Pixel 57	18	Pixel 58	33	AGnd / KGnd
4	Pixel 55	19	Pixel 56	34	AGnd / KGnd
5	Pixel 53	20	Pixel 54	35	AGnd / KGnd
6	Pixel 51	21	Pixel 52	36	AGnd / KGnd
7	Pixel 49	22	Pixel 50	37	AGnd / KGnd
8	Pixel 47	23	Pixel 48	38	AGnd / KGnd
9	Pixel 45	24	Pixel 46	39	AGnd / KGnd
10	Pixel 43	25	Pixel 44	40	AGnd / KGnd
11	Pixel 41	26	Pixel 42	41	AGnd / KGnd
12	Pixel 39	27	Pixel 40	42	AGnd / KGnd
13	Pixel 37	28	Pixel 38	43	n/c
14	Pixel 35	29	Pixel 36	44	Pixel 34
15	Pixel 33	30	Pixel 32	-	-

## HD44 #3

1	Pixel 92	16	n/c	31	Test pixel 3
2	Pixel 91	17	Pixel 93	32	n/c
3	Pixel 89	18	Pixel 90	33	AGnd / KGnd
4	Pixel 87	19	Pixel 88	34	AGnd / KGnd
5	Pixel 85	20	Pixel 86	35	AGnd / KGnd
6	Pixel 83	21	Pixel 84	36	AGnd / KGnd
7	Pixel 81	22	Pixel 82	37	AGnd / KGnd
8	Pixel 79	23	Pixel 80	38	AGnd / KGnd
9	Pixel 77	24	Pixel 78	39	AGnd / KGnd
10	Pixel 75	25	Pixel 76	40	AGnd / KGnd
11	Pixel 73	26	Pixel 74	41	AGnd / KGnd
12	Pixel 71	27	Pixel 72	42	AGnd / KGnd
13	Pixel 69	28	Pixel 70	43	n/c
14	Pixel 67	29	Pixel 68	44	Pixel 66
15	Pixel 65	30	Pixel 64	-	-

**PSI System Controls and Diagnostics** 



#### Pixel readout (cont) HD44 #4. 1 Pixel 124 16 n/c Test pixel 4 31 Pixel 125 2 Pixel 123 17 32 n/c 3 Pixel 121 18 Pixel 122 33 AGnd / KGnd AGnd / KGnd 4 Pixel 119 19 Pixel 120 34 5 Pixel 117 20 Pixel 118 35 AGnd / KGnd 21 Pixel 116 AGnd / KGnd 6 Pixel 115 36 7 Pixel 113 22 Pixel 114 37 AGnd / KGnd Pixel 112 AGnd / KGnd 8 Pixel 111 23 38 Pixel 109 Pixel 110 AGnd / KGnd 9 24 39 10 Pixel 107 25 Pixel 108 40 AGnd / KGnd 11 Pixel 105 26 Pixel 106 41 AGnd / KGnd 12 Pixel 103 27 Pixel 104 42 AGnd / KGnd

HV in	SHV
HV out	SHV
	HV in and out are interchangeable.

28

29

30

Pixel 102

Pixel 100

Pixel 96

43

44

n/c

Pixel 98

13

14

15


Pixel 101


Pixel 99


Pixel 97

Connectors (cont)

## Calibration Gain curves Approximate gain curves at standard ambient temperature and pressure for protons, 7.62 mm gap. 400.0 350.0 300.0 250.0 200.0 150.0 100.0 50.0 0.0 300 Proton energy (MeV) 7.62 mm 200.0 150.0 Gain 100.0 50.0 100 150 200 250 Proton energy (MeV) \_\_\_\_\_7.62 mm Note: Critical dosimetry measurements must use accurate gain values referenced to traceable standards, and regularly validated.







Side elevations Dims mm

## **Ordering information**

PX-2

Pixelated ionization chamber with 4.2 cm diameter sensitive area, thin FR4 cathode with 120 pixels.

Pyramid Technical Consultants, Inc., 1050 Waltham Street Suite 200 Lexington MA 02421 USA Tel: +1 781 402 1700 (USA),

+44 1273 492001(UK)

Email: support@ptcusa.com

www.ptcusa.com

The information herein is believed accurate at time of publication, but no specific warranty is given regarding its use. All specifications are subject to change.

All trademarks and names acknowledged.

PX-2\_DS\_150507

**PSI System Controls and Diagnostics** 



**Pyramid Technical Consultants**